

## Voluntary Recall of Selected Lots of ACUVUE® Brand Contact Lenses

At Johnson & Johnson Vision, our top priority is patient safety and we hold ourselves to high standards for product quality and customer satisfaction.

We are voluntarily recalling limited lots of ACUVUE® Brand Contact Lenses in the United States and Canada because they do not meet quality standards.

Importantly, no serious adverse events have been reported due to these issues. We have already taken measures to help ensure these issues do not recur. ACUVUE® Brand Contact Lenses are safe when used as directed and can continue to be used with confidence.

You can find affected product lots and names in the chart below, as well as instructions for finding the lot information on your ACUVUE® Brand Contact Lenses.

Reasons noted do not affect all lenses in these lots. However, our recall is for the full lots listed below.

Brand Name	Lot Number	Parameters	Reason for Recall
ACUVUE® ADVANCE®	<b>B00DHLP</b>	BC: 8.7 Power: -10.00	Potential for brush bristle to be found between the blister package and foil
ACUVUE® 1-Day MOIST for ASTIGMATISM	<b>B00LF41</b> (30-pack) <b>F00LF41</b> (90-pack)	BC: 8.5 Power: -7.50 Cyl: -1.25 Axis: 180	
ACUVUE OASYS® for ASTIGMATISM	<b>B00GW4Z</b> (6-pack) <b>B00GW4ZC</b> (1-pack)	BC: 8.6 Power: -6.00 Cyl: -1.75 Axis: 170	Lens may not provide appropriate vision correction "power," which may result in some distortion or blurriness.
ACUVUE OASYS® for ASTIGMATISM	<b>B00HRMG</b> (6-pack) <b>B00HRMG4</b> (3-pack)	BC: 8.6 Power: -3.50 Cyl: -0.75 Axis: 180	

The lot number is displayed in the barcode area on the back of each individual unit carton (below for example). The lot number is also present on the foil of each individual blister package of the contact lens.

Johnson & Johnson Vision Care, Inc. has voluntarily initiated this recall. We have identified the root causes and confirmed that improvements at our manufacturing sites have been completed as corrective actions.

If you have a product that might be affected, please **take the following actions**:


1. **Review** your inventory and determine if you have ACUVUE® Brand Contact Lenses **with any of the Lot Numbers listed above.**
2. **STOP** using all **affected** product. Note: You can continue to use product from all other lots not affected by this recall.
3. **Contact** Johnson & Johnson Vision Customer Relations at 1-800-843-2020 to arrange for return of your product and to provide you with replacement product.

If you have product complaints or adverse events to report regarding the use of ACUVUE® Brand Contact Lenses, please call 1-800-843-2020. If you do report a complaint, please provide the ACUVUE® Brand Contact Lenses lot number and a description of the issue.

\*\*\*\*\*

**Product Unit Carton & Blister Pack Label Examples**  
**(Actual packaging appearance may vary by product and lens specifications)**

(ACUVUE OASYS® for ASTIGMATISM)



Example: Lot Number B00HJ3H

**(Carton)**

Barcode area, right side of **LOT** of reverse side of the carton, seven digits of alphanumeric from the left.



**(Blister Package)**

Bottom side of the foil, right side of **LOT**, seven digits of alphanumeric.

